



SECHRIST HYDRAULIC GURNEY LOW PROFILE MODEL


PRODUCT FEATURES:

Sechrist has the unique distinction of offering the only Low Profile Gurney currently in the market. Compatible with all Sechrist Monoplace Hyperbaric Chambers, this new product offers a variety of features and benefits that will meet the needs of your Hyperbaric patients and Hyperbaric team.

The Sechrist Low Profile Gurney is provided with the H Series Monoplace Hyperbaric Chamber Storage Models. Convenient storage allows for additional free space in the hyperbaric room eliminating need for wall space to store chamber gurneys.

28" from mattress to floor

23" from mattress
to floor



SECHRIST LOW PROFILE GURNEY



*Low Profile Gurney (on left) compared to
Standard Gurney (on right)*

- Low Profile Gurney goes 5" lower from floor to top of mattress than Standard Gurney allowing for maximum patient comfort and safety.
- Eliminates or minimizes the use of step stools and staff assistance which potentially results in easier and safer patient access to getting onto the Low Profile Gurney.
- Adjustable backrest for Models 3600H and 4100H has an increased incline up to 45° for additional flexibility and patient comfort.
- Bariatric patients up to 700 lbs can be comfortably supported.
- The Low Profile Gurney provides for easy loading of the patient into the chamber without additional transfers.

Phone: 800-SECHRIST (800-732-4747) • 714-579-8400

Corporate Headquarters: 4225 E. La Palma Ave. • Anaheim • CA 92807

Website Addresses: www.sechristusa.com.com | www.sechriststretchers.com

- Sechrist Low Profile Gurneys offer considerable flexibility, easily adjusting to various heights for ready access to patients in beds, wheelchairs or other elevations.
- Designed exclusively for all Sechrist Monoplace Hyperbaric Chambers allowing for easy transferring of patients in and out of Sechrist Hyperbaric Chambers.
- One step locking mechanism ensures simultaneously locking of all wheels for safety.
- Added features include side rails for patient protection, IV pole mount, and storage tray.


AUTOMATIC GURNEY MECHANICAL HEIGHT STOP

Sechrist's exclusive Automatic Gurney Mechanical Height Stop will allow users to always stop the gurney to meet the chamber rails at exactly the right height every time.

PRODUCT BENEFITS:

- Efficiency and simplification of patient treatments not only enhances the patient experience but also improves the process for the hyperbaric technician
- Reduces the need to continually check to see if the rails are aligned
- No longer requires 'feathering' of the gurney height which can lead to stretcher plate becoming uneven
- Difficulty in lining up gurney to chamber is eliminated
- Minimizes potential damage to chamber and door gasket when aligning

All new Low Profile Gurneys come with the Automatic Gurney Mechanical Height Stop. Available as an upgrade to most older Sechrist Low Profile Gurneys, a retrofit kit can be conveniently installed at the time of Annual Service.



ORDERING INFORMATION

PRODUCT	PART NUMBER
OmniPressure Mattress	20505
OmniPressure Bariatric Mattress	20509
Track Guide Bumpers (10 pack)	HB601-06D
Stretcher Assembly	21710
Stretcher Plate Conversion Kit	20562
Automatic Gurney Mechanical Height Stop Retrofit Kit	20655
Low Profile Gurney Base	21709-05

P/N 100459, rev 6 (March 2024)



Registered Address:
Obelis s.a.
Bd. Général Wahis 53
B-1030 Brussels, Belgium
E-mail: mail@obelis.net