

SECHRIST
CERTIFIED
PRE-OWNED
(CPO)
MTI CHAIRS



SECHRIST CERTIFIED PRE-OWNED (CPO) MTI CHAIRS

MTI MODELS 527,
527W, 529 AND 529W

RESTORED MTI CHAIRS



MAIN FEATURES:

- **CONTROLS:** Foot control with optional foot section mounted hand control
- **PROGRAMMABLE BUTTONS:** 4 dynamic user-programmable buttons
- **HOME BUTTON:** ADA compliant home and user-programmable home positions
- **HEADREST:** Integrated into backrest
- **PATIENT ARMS:** Slide-back, floating and removable
- **LEVELING GLIDES:** Adjustable and built-in
- **G2 ACCESSORIES:** Plug and play accessories

WHY MTI CHAIRS?

MTI's Americans with Disabilities Act (ADA) compliant accessible chairs can help reduce patient and healthcare worker injury and distress to achieve a true patient and healthcare worker centered design. They also help reduce liability issues and increase patient satisfaction and employee retention using the newest technology and quality standards. MTI Models 529 and 529W only comply with 36 CFR Part 1195 and 2010 ADA Standard when used with optional transfer support. *(Source MTI).*

GET IN TOUCH TODAY
FOR MORE INFO



www.sechristusa.com
www.sechristcpo.com



P : 714-579-8400
E : info@sechristusa.com

4225 E. La Palma Ave
Anaheim, CA, 92807

MODEL COMPARISON

MODEL	527	527W
Power Functions	Power lift, back and tilt functions (non-knee break)	Power lift, back and tilt functions (non-knee break)
Programming	Standard two (2) position programmability	Standard two (2) position programmability
Patient Capacity	650lb	800 lb
Seat Width	24"	30"
Seat Entry - Extension Height	19" entry height with 20" vertical travel power lift base	19" entry height with 20" vertical travel power lift base
Base	Standard	Standard
Standard Accessories	Foot control with auto return, stop/cancel, and any button cancellation	Foot control with auto return, stop/cancel, and any button cancellation
ADA Compliant	Yes	Yes

MODEL	529	529W
Power Functions	4 - Lift, Back, Tilt, Foot Ext.	4 - Lift, Back, Tilt, Foot Ext.
Programming	4 Program Buttons, Programmable Home & ADA Compliant Home	4 Program Buttons, Programmable Home & ADA Compliant Home
Patient Capacity	725 lb	850 lb
Seat Width	24"	30"
Seat Entry - Extension Height	18.75-39" on Standard Base 18.75-37" on Swivel Base	18.75-39"
Base	Standard or Swivel	Standard
Standard Accessories	Programmable Foot Control and Slide-back Patient Arms	Programmable Foot Control and Slide-back Patient Arms
ADA Compliant	Yes	Yes

SMART TECH™ FEATURES

SmartTech™ is the combination of MTI's Smart Controller™ and Smart Safety™ technology using a centralized computer controller that enhances safety, reduces energy consumption, and provides smoother movement, additional memory options and more user preferences.



ADA

Compliant with All ADA Regulations¹



CONTROL

Foot Control with Sealed Membrane for Cleaning and Disinfecting



EXTENSION

Power Footrest Extends 11" for Taller Patients



HASSLE FREE

Maintenance-Free Joints. No Lubricating Pivot Joints... EVER!



HIBERNATE

After 60 Minutes Chair Sleeps to Save Energy



SPEED

Auto-Lock and Lock/Unlock Electronic Swivel Settings (Optional)



SWIVEL

Maintenance-Free Joints. No Lubricating Pivot Joints... EVER!



HOME

ADA Compliant Home and Programmable Home Positions



MEMORY

4 Dynamic User-Programmable Buttons



MOVEMENT

Removable Slide Back Patient Arms Compensate with Backrest Movement



SAFETY

Smart Safety™ Technology Helps Keep the Chair in Safe Positions



SOFT START

Technology Providing Smooth Start and Stop Movements



UNIVERSAL

Universal Power Supply for 100-240 Volt Applications (VAC)



UPHOLSTERY

Quick Installation and Removable for Easy Cleaning and Replacement

¹36 CFR Part 1195 and 2010 ADA Standard when used with optional transfer support



DISCOVER THE ADVANTAGES OF SECHRIST CERTIFIED PRE-OWNED (CPO) MTI CHAIRS



WHY SECHRIST PRE-OWNED?

The work we put into refurbishing Certified Pre-Owned (CPO) MTI Chairs makes us the perfect choice for those seeking exceptional value without compromising on the renowned quality, innovation, safety, reliability, and service associated with Sechrist.

QUALITY COMMITMENT:

When you opt for Sechrist Certified Pre-Owned MTI Chairs, you're investing in products that have undergone meticulous inspection and servicing by our skilled technicians. We perform a thorough assessment, guaranteeing you a chair that adheres to the most stringent quality benchmarks.

GENUINE PARTS

We use only genuine, authorized MTI replacement parts that adhere to the strict quality standards set by Sechrist Industries.

COMPREHENSIVE SUPPORT:

When you choose a Sechrist Certified Pre-Owned MTI Chair, you'll have access to a wide spectrum of services, including technical guidance, and ongoing support for your chair from Sechrist technicians trained and certified by MTI.

WARRANTY AND MAINTENANCE:

All Sechrist Certified Pre-Owned MTI Chairs come with a 90-day limited warranty. Plus, we offer a range of annual service and maintenance program options. Our dedicated technical team is available 24/7 to assist you.

Unlock the opportunity to own a Sechrist CPO MTI chair for less than you might have imagined. Discover the advantages of Certified Pre-Owned today! Logistics, delivery, installation, repair and annual maintenance performed by Sechrist. We look forward to assisting you.

www.sechristusa.com
www.sechristcpo.com



P : 714-579-8400
E : info@sechristusa.com

4225 E. La Palma Ave
Anaheim, CA, 92807