



INDUSTRIES, INC.



# SECHRIST WALL MOUNT GROUND KIT

ENSURING THE PROPER PREPARATION OF THE HYPERBARIC SUITE IS KEY TO THE SUCCESSFUL INSTALLATION OF SECHRIST MONOPLACE HYPERBARIC CHAMBERS. GROUNDING OF SECHRIST'S MONOPLACE HYPERBARIC CHAMBERS IS A VITAL SAFETY CONSIDERATION AND NFPA HEALTHCARE FACILITIES CODE 99 REQUIREMENT FOR CLASS B CHAMBERS.

## PRODUCT FEATURES:

Kit is complete with everything needed to properly install the Wall Mount Ground for Sechrist Hyperbaric Chambers

- Wall Plate
- Threaded Ground Stud
- Hardware and Installation Information

**Wall Mount Ground Kit meets all Sechrist specifications per the Technical Installation Guide. Ordering Information - P/N 20656**

## HEADQUARTERS

4225 E. La Palma Ave. - Anaheim, CA 92807  
Phone : 1-800-SECHRIST (800-732-4747) • Fax : 714-579-0814  
customerservice@sechristusa.com • www.sechristusa.com

P/N 100582, REV 0 (DEC 2022)