




INDUSTRIES, INC.

ENTERTAINMENT SYSTEM MOUNTING KIT MEDICAL GRADE & STANDARD TV SYSTEMS

*TV not included with kit;
for display purposes only



Many of our hospital and clinic partners would like to utilize their own TV with the Sechrist Entertainment System. This option is available and made convenient through the Sechrist Entertainment Mounting Kit which includes all the necessary parts and equipment to accommodate both Medical Grade TV Systems and Standard TV Systems (power cord sold separately). Easy to follow instructions allow for simple installation by either a Sechrist Authorized Technician or hospital/clinic biomed personnel. The Mounting Kit provides flexibility when selecting the Entertainment System for your hyperbaric patients.



*TV not included with kit; for display purposes only

THE WORLD LEADER IN HYPERBARIC TECHNOLOGY & MEDICINE

ENTERTAINMENT SYSTEM MOUNTING KIT

MEDICAL GRADE & STANDARD TV SYSTEMS

| MOUNTING KITS | | PART NUMBER |
|--|---|-------------|
| FOR STANDARD TV SYSTEM | | |
| Models 3200 / 3300E / 3300H / 3600E / 3600H Series Hyperbaric Chambers | | 21868 |
| Model 4100H Series Hyperbaric Chambers | | 21866 |
| FOR MEDICAL GRADE TV SYSTEM | | |
| Models 3200 / 3300E / 3300H / 3600E / 3600H Series Hyperbaric Chambers | | 21869 |
| Model 4100H Series Hyperbaric Chambers | | 21867 |
| TECHNICAL RECOMMENDATIONS PER SYSTEM | | |
| TV SIZE | TV must not exceed 24 inches | |
| MOUNTING PATTERN | TV mounting pattern must be Versa universal mount 100mm x100mm; usually standard for 24 in and small TVs | |
| AUDIO OUTPUT STANDARD TV | Audio output must have 3.5mm or 1/4 in. headphone jack output if utilizing Sechrist audio output cable | |
| AUDIO OUTPUT MEDICAL GRADE TV | Audio output must have 6-pin pillow port output if utilizing Sechrist audio output cable | |
| POWER CORD | TV power cord must be sufficiently long to reach from outlet to top of chamber**, 12 ft. recommended | |
| KIT COMPONENTS AND PART NUMBERS | | |
| AUDIO OUTPUT CABLE FOR STANDARD TV | TV audio output standard cable 3.5mm headphone jack to male dual RCA; (P/N 22277) Qty 1, 15ft | |
| AUDIO OUTPUT CABLE FOR MEDICAL GRADE TV | TV audio output cable 6-pin pillow port to male single RCA; (P/N 22248) Qty 1, 12 ft. | |
| AUDIO CABLE | Audio cable: 1/4 in. phone plug to male dual RCA; (P/N 22247) Qty 1, 12 ft. | |
| CABLES & HARDWARE | Coaxial Cable, Video low profile; (P/N 22245) Qty 1, 25 ft. Tubing Split Corrugated 3/4 in. - inner diameter (P/N 43541) Qty 1, 12 ft. Cable Holder Adhesive back (P/N 43715) Qty 5 Ground Loop Isolator (P/N 41063) Qty 1 Screw, Panhead Phillips M4 x 6 mm Long (P/N 56541) and 8mm Long (P/N 46541) Qty 4 each | |
| POWER CORD (not included in kit) ** | Power Cord, 13A, 125VAC, IEC320 female receptacle; (P/N 22223), 15 ft **Note - please verify that cord receptacle fits TV purchased (see image below). | |



****Shown: Power Cord Receptacle (Part number 22223)**

sechristusa.com

THE WORLD LEADER IN HYPERBARIC TECHNOLOGY & MEDICINE